

Excellence Weighing Solutions



Weighing Solutions Guide

Innovative features

Built-in applications

Useful accessories

Weighing Solutions

For more Security and Productivity

METTLER TOLEDO

Weighing Solutions Guide for METTLER TOLEDO Excellence balances

The Excellence balance range from METTLER TOLEDO offers many innovative features, built-in applications and useful accessories which have been carefully designed to meet the needs of our customers around the world. One of which should certainly provide the right answer to your specific weighing process. This "Solution Guide" introduces some of our unique and innovative solutions - our answers to your daily challenges!



Table of Contents

Features overview of the Excellence balance range _____ 4

Your Challenges...

Our Solutions...

1. Regulatory Compliance

5

1.1. Does your SOP require you to level the balance at regular intervals?	LevelControl _____	6
1.2. Do you have to fulfil USP minimum weight requirements?	MinWeigh _____	7
1.3. Do you need different user authorization levels?	User Management _____	8
1.4. Do you have to comply with 21 CFR part 11?	LabX balance pro _____	9

2. More Security

11

2.1. Are you troubled by static electricity?	AntiStatic kit _____	12
2.2. Do you want to weigh directly into various tare containers?	ErgoClips _____	14
2.3. Do you ever weigh hazardous materials?	SmartSens & ErgoSens _____	15
2.4. Do you weigh in a hazardous environment?	Ex balances _____	16
2.5. Is your balance always adjusted correctly?	proFACT / FACT & BalanceCheck _____	17
2.6. Do you have cross-contamination risks?	Cleaning concept for Excellence balances _____	18

3. Increased Productivity

19

3.1. Would you like to increase the speed of your weighings?	SmartGrid & Draft shield for XP/XS Analytical balances _____	20
3.2. Would you like to improve the ease of use of the balance?	Touch screen operation _____	21
3.3. Are you looking for software to manage weighing data?	LabX balance _____	22
3.4. Do you ever weigh magnetic samples?	MPS weighing pan _____	23
3.5. Are cables making a mess of your laboratory?	Bluetooth wireless connection _____	24
3.6. Do you need a customized draft shield?	Draft shields for XP/XS Precision balances _____	25
3.7. Future software upgrades?	e-Loader _____	26
3.8. How soon can you expect to make a return on your investment?	ROI calculation example _____	27

4. More Flexible Connectivity

29

4.1. Do you connect your balance with other peripherals?	Optional interfaces for Excellence series _____	30
	Connectivity with peripherals _____	31
	MT-SICS commands _____	32
4.2. Are you about to replace your old balance?	Replacement models _____	33
	Backwards compatibility _____	33
4.3. Do you want to use a USB interface?	USB connections _____	34
4.4. Do you connect your balance to a titrator?	Titrator connection _____	35

5. Comprehensive Range of Accessories

37

5.1. Accessories overview for Excellence series _____	38
---	----

Excellence Plus XP


Features overview of the Excellence balance range

XP Analytical balances

- MinWeigh*
- proFACT
- User management (8 users)
- GxP compliance
- SmartScreen (color touch screen display)
- QM tool box
- LevelControl
- BalanceCheck
- SmartSens
- SmartGrid
- AntiStatic kit*
- ErgoClips*
- Connectivity*
- LabX balance - Data & instrument management*

Built-in applications

- Weighing
- Statistics
- Formulation
- Piece counting
- Percent weighing
- Density determination
- Differential weighing




XP Precision balances

- MinWeigh*
- proFACT
- User management (8 users)
- GxP compliance
- SmartScreen (color touch screen display)
- QM tool box
- LevelControl
- BalanceCheck
- SmartSens
- MPS weighing pan*
- Connectivity*
- LabX balance - Data & instrument management*

Built-in applications

- Weighing
- Statistics
- Formulation
- Piece counting
- Percent weighing
- Density determination
- Differential weighing
- Dynamic weighing




Excellence XS

XS Analytical balances

- FACT
- GxP compliance
- Touch screen
- SmartGrid
- ErgoClips*
- Connectivity*
- LabX balance - Data & instrument management*

Built-in applications

- Weighing
- Statistics
- Formulation
- Piece counting
- Percent weighing
- Density determination




XS Precision balances


- FACT
- GxP compliance
- Touch screen
- MPS weighing pan*
- Connectivity*
- LabX balance - Data & instrument management*

Built-in applications

- Weighing
- Statistics
- Formulation
- Piece counting
- Percent weighing
- Density determination
- Dynamic weighing



* Optional accessories / services



1

1. Regulatory Compliance

1.1 Does your SOP require you to level the balance at regular intervals?

Could you imagine a balance that has a self-checking system for level adjustment?
Our **Excellence Plus XP** is equipped with a **LevelControl** system, which automatically warns when the balance is not leveled correctly.

Your challenges


- A balance can be moved out of level for many reasons.
- If the balance is not correctly leveled, the performance of the balance cannot be guaranteed (sensitivity offset is increased).
- Companies working to guidelines such as ISO and USP may be required to check that the balance is leveled each time before use.




An optical sensor is always keeping eye on the level.

Our solutions

- The **LevelControl** system of the **Excellence Plus XP** automatically detects if the balance is leveled or not.
- When the balance is inclined, **LevelControl** gives a message and/or an acoustic warning.
- The balance gives you leveling instructions on the display.
- The leveling status is always stored in the Adjustment History.
- The leveling status can be printed in the header, in the results or with each single weight value, in every application.



When the balance is out of level, a warning message appears.



An onscreen graphic shows you how to rapidly level your balance.

```

----- Weighing -----
10.May 2006      18:54
User Name      Home
Balance Type    XP205DR
WeighBridge SNR: 12345678
Terminal SNR: 1234567890

Balance is levelled
                  0.00000 g

Balance is levelled
                  7.58942 g

Balance is inclined
                  7.58946 g
  
```

Printout example

Your benefits

- In the daily weighing process, you no longer need to check whether the balance is leveled.
- The leveling status can be printed out or stored in the balance memory. The traceability of instruments and the validity of weighing results are guaranteed.

Availability

The **LevelControl** feature is standard on all **XP Analytical** and **XP Precision balances**.

1.2 Do you have to fulfil USP* minimum weight requirements?

Your challenges


- When you dispense **small amounts of sample**, are you sure your weighing results are always correct? Sometimes weighing processes are carried out at between 1 % and 10 % of a balance's weighing capacity. But it is precisely in this range that the greatest relative measurement uncertainties occur.
- If the weighing value is less than the **minimum weight**, the weighing result cannot be guaranteed to be within process requirements. In a regulated environment, the measurement series must be repeated, and valuable substances may be lost. In the worst case, invalid values could cause production stops, costing millions.
- **Minimum weight** depends on the requirement of each regulation (e.g. USP regulation requires an accuracy of better than 0.1 % with 3 times the standard deviation).

Our solutions


- The **MinWeigh application** of the **Excellence Plus XP** determines the **minimum weight** according to your compliance regulations.
- **Minimum weight** is determined in the actual working environment where the balance is installed.
- If the value is below the determined **minimum weight**, the color of the displayed value is **red** to warn that the current value is not sufficient to ensure tolerances are met.
- **An asterisk** is printed out next to the result if the value is below the **minimum weight**.
- Up to **3 different methods** can be determined and stored. Every user can therefore recall the unique method which is required.

1.2

MinWeigh



When the displayed value is below minimum weight, it is **red** to signal a violation.



Up to 3 different methods can be configured.

----- Weighing -----	
User Name	Home
MW-Method	GLP
RefTare	1.00000 g
MinWeigh	0.01000 g
T	1.01787 g
N	1.00145 g
G	2.01932 g

Printout example

Your benefits

- With **MinWeigh** you can always be sure that each result fulfills the requirements, preventing out of spec results.
- **MinWeigh** saves time and prevents waste of precious substances. Saving money, and ensuring compliance.

Availability

MinWeigh is a built-in application for **XP Analytical** and **XP Precision balances**. It must be activated by a certified MT service technician at the installation site. After programming the balance, the service technician issues a "MinWeigh Certificate" documenting the test measurements and tolerances, as well as the corresponding tare and minimum weight values.

*USP: United States Pharmacopeia

1.3 Do you need different user authorization levels?

Your challenges

- Weighing results are often just the first stage of a much larger laboratory process. For regulatory compliance, changes to an instrument's settings might influence the validity of the weighing results and so must be strictly controlled. And each change to an instrument must be documented for traceability of the instruments change history.
- In a regulated laboratory, the administrative rights are described for each user in detail.




Depending on user level, detailed individual access rights can be configured (by the administrator).



Traceability is guaranteed with the Change History.

Our solutions

- Individual user settings can be configured for **up to 8 users**. The "Administrator" can switch off any unnecessary user settings.
- The **User Management** function of the **Excellence Plus XP** enables the "Administrator" to control each user's access rights. Those access rights can be flexibly configured depending on each user's level. (E.g. User A can execute an adjustment using an external weight or can use different applications like "Weighing" or "Statistics" freely. Whereas User B can only use one application and is not permitted to make any changes to parameters or instrument settings.)
- Individual access rights can be set by the administrator, depending on the responsibilities of the users.
- Each user setting can be protected using a password.
- The "Administrator" can also switch off unnecessary applications in each user setting.
- In the change history, the last 50 actions are automatically stored, and can be printed out as required.



Up to 8 users are available in the system and unnecessary users can be switched off by the "Administrator".

Your benefits

- With the **User Management** function, you can rule out any unauthorized operation or changes to equipment settings.
- Traceability of the balance is guaranteed.

Availability

The **User Management** function is available in all **XP Analytical** and **Precision balances**.

XS analytical and precision balances have only 1 user. Simple password protection of the balance is possible.

1.4 Do you have to comply to 21 CFR part 11?

Your challenges

- You need laboratory software to connect balances in your laboratory to a network. The software must be connected with a LIMS for comprehensive data management.
- Instruments need to be controlled remotely by the software, and all data needs to be checked directly from the office. At the same time, **audit trails, electronic signatures and approvals/sign-off of documents** must be fulfilled for **21 CFR part 11 compliance**.



LabX pro guarantees complete data management in the laboratory and supports **21 CFR part 11** requirements.

Our solutions & your benefits

- With **LabX pro** there is no need to go to the laboratory to check results. Just start LabX pro on your office PC and you can check results, change methods and review and sign off documents electronically; you don't have to print a single document. And it supports **21 CFR part 11**.
- With **LabX pro**, you can simplify your tasks in many ways: for example you can import your sample data directly from the LIMS, make templates for routine daily tasks, create statistics and control charts automatically.
- Using Rapid Access in **LabX** you do not need to touch a computer. **Excellence balances** allow you to log on to LabX directly from the balance. There is no need to connect the balance to a computer either: Just plug it directly into your company network with the optional **Ethernet interface**.


1.4

LabX pro

Availability and functionality overview

Rapid Access with Excellence Plus XP balances - as if there were no computer.

1



With the **Rapid Access buttons**, select the job you want to carry out from the job queue.

2



Log on to LabX from the balance touch screen. The balance prompts you for your user ID and password.

3



LabX checks the release status of the balance. You can then follow the **step-by-step** guidance through the job on the balance display.

4



Press the print button to transfer data.

5



Retrieve data, get a **printed report** or have data **automatically exported** to a LIMS with LabX connect.

LabX is the laboratory instrument control software from METTLER TOLEDO, for **Balances** and **Titration**.

There are 3 levels: **LabX pro**, **light** and **direct**.

	21 CFR Part 11 support	LIMS integration	Rapid Access*	Data management	Remarks
LabX pro balance	X	X	X	X	Full compliance
LabX pro titration	X	X	n/a	X	Full compliance
LabX light balance	n/a	n/a	X	X	Complete data management
LabX light titration	n/a	n/a	n/a	X	Complete data management
LabX direct balance	n/a	n/a	n/a	n/a	Direct data transfer to Excel, Word etc.

LabX balance can control METTLER TOLEDO XP/XS, AX/MX/UMX, AT/MT/UMT, AB-S/(FACT), PB-S/(FACT), AG, PG-S, SG, PR, SR balances.

LabX titration can control METTLER TOLEDO DL50 Graphix/Rondolino, DL53, DL55, DL58, DL70ES, DL77, DL31, DL32, DL38, DL39 titration systems.

For more detailed information

Please visit www.mt.com/LabX.

[illegible]



2

2. More Security

2.1 Are you troubled by static electricity?

Static electricity is one of the major factors affecting good weighing results.

Are you weighing different powders using plastic or glass containers?

You might have to consider the effects of static electricity.

The optional AntiStatic kit for the XP Analytical offers the most effective and fastest solution to ensure **accurate weighing results**, even from critical or lightly charged samples.

Your challenges

- Insulating materials such as plastic or glass are affected easily by static electricity, which can generate an **electrostatic force**. This force causes unstable weighing results; the display value may keep changing and never settles.
- The weighing result is not only unstable, but might also be wrong, since static forces can cause weighing errors of plus or minus several milligrams!
- Powdery samples tend to be scattered or stick to the tare container due to static electricity. It is not only difficult to handle, but it can also lead to cross-contamination. If the sample is toxic or harmful, this could be dangerous to the operator.


Positive & negative ions are generated, **without air turbulence!**

Our solutions


- For **XP Analytical balances**, the **built-in AntiStatic kit** is available as an optional accessory.
- The **AntiStatic kit** can be fitted in six different positions to the draft shield (depending on the sample size and shape, and whether right or left handed).
- Ionization can be activated automatically in response to opening the draft shield door or can be set manually to work as required.

2.1

AntiStatic kit



The **AntiStatic kit** can be positioned at 3 different heights*, **left side or right side** (6 places in total).
The middle position is recommended for most cases.



Ionization is synchronized with the draft shield movement. It is activated when the door is opened and/or after the door is closed. Both ways can be configured flexibly, together with time activation (we recommend to set the "ionization time" less than 5 seconds for faster weighing result).

Your benefits

- The **AntiStatic kit** for **XP Analytical** generates positive and negative ions evenly, and due to its high voltage the static electricity of the sample or tare containers can be discharged **within a few seconds**.
- The **AntiStatic kit** **does not** generate air turbulence, which could cause powdery samples to rise into the air. The risk of cross-contamination is therefore reduced to a minimum.
- Thanks to complete integration with the balance, the antistatic kit is activated automatically with door movement (i.e. you don't have to operate the **AntiStatic kit** manually for every weighing).
- This results in extremely quick, efficient and even more accurate weighings, **with maximum security of sample, user and results**.

Availability and ordering information

The **built-in AntiStatic kit** can only be installed with the **XP Analytical balance** as an optional accessory.

Accessories	Description	Part No.
AntiStatic kit	Built-in AntiStatic kit for XP Analytical balance, with power supply	11107761
	Optional second point-electrode for AntiStatic kit	11107762

Tips and recommendations

- For the Excellence XS Analytical balance or other balances, we recommend using the **Universal AntiStatic kit**. This also generates positive and negative ions evenly, and causes no air turbulence.




Universal AntiStatic kit can be placed freely.

Accessories	Description	Part No.
Universal AntiStatic kit	Universal AntiStatic kit complete (U-shaped), including power supply	11107767
Optional U-shaped electrode	U-electrode for Universal AntiStatic kit	11107764
Optional point electrode	Point electrode for discharging small samples	11107765
Power supply	Power supply for Universal AntiStatic kit (required for U-shaped electrode and point electrode)	11107766

2.1

AntiStatic kit

- To reduce the influence of static electricity, covering a statically charged sample or tare container with a metal cage is also effective, due to the "**Faraday effect**".
- Our **ErgoClips** offer an additional help for stable results. The **AntiStatic kit + ErgoClips** offer the best solution against static electricity.
- For more detailed information on **ErgoClips**, please refer to chapter 2.2.



The electrostatic force lines are screened by the basket.



ErgoClips reduce the influence of static electricity.

2.2 Do you want to weigh directly into various tare containers?

Your challenges

- Some typical laboratory tare containers are difficult for weighing because of their odd shape.
- Samples are sometimes spilled because of unstable and/or narrow tare containers.
- One step dosing directly into the tare container offers less contamination risk and minimizes loss of sample, enabling faster dosing procedures.

Our solutions & your benefits

- We offer various type of weighing accessories called **"ErgoClips"**, which can be fixed to the grid weighing pan on **XP/XS Analytical balances** in order to make your weighings faster and easier.
- **ErgoClips** hold tare containers securely and allow you one-step dosing, directly into the tare container on the balance.
- The **ErgoClips** can reduce the influence of static electricity, thanks to the **"Faraday effect"**.

Availability and ordering information

ErgoClips are fitted to the grid pan of **XP/XS Analytical balances**.



ErgoClip Basket small
Part No. 11106747

1

Suitable for several test tubes and syringes of different sizes. The basket shape reduces the effects of static electricity. **This ErgoClip is included with all XP/XS Analytical balances.**



ErgoClip Titration Basket
Part No. 11106746

2

Suitable for bigger tare containers, like Titration beakers. The large basket reduces static electricity effects and ensures secure placement.



ErgoClip Weighing Boat
Part No. 11106748

3

Suitable for lab weighing boats. The basket shape reduces static electricity effects.



ErgoClip Round Bottom Flask
Part No. 11106746

4

Suitable for various sizes of round bottom flasks, to guarantee secure placement.



ErgoClip Flask
Order No. 11106764

5

Specialized ErgoClip for graduated flasks. Due to less surface area, the weighing accuracy is improved, and stabilization time is very short. Suitable for dosing small amount of sample directly into a flask.



MinWeigh Door
Part No. 11106749

6

Optional draft shield with a small adjustable opening. Perfect combination with the ErgoClip Flask (5), for even faster stabilization time. The MinWeigh Door allows the minimum weight to be reduced even further.



ErgoClip Flask micro
Part No. 11107879

7

Specialized ErgoClip for small graduated flasks. Due to less surface area, the weighing accuracy is improved, and stabilization time is very short. Suitable for dosing small amounts of sample directly into a flask (for XP micro balances only).



MinWeigh Door micro
Part No. 11107869

8

Optional inner draft shield with a small adjustable opening. Perfect combination with the ErgoClip Flask micro (7), for even faster stabilization time. The MinWeigh Door micro allows the minimum weight to be reduced even further (for XP micro balances only).



ErgoClip Tube
Part No. 11106784

9

Specialized ErgoClip for test tubes. Due to less surface area, the stabilization time is very short. Suitable for dosing small amounts of sample directly into the tube.



Weighing Kit - Set of ErgoClips
Part No. 11106707

10

Set of ErgoClips 1, 3, 4, disposable aluminum pans, and weighing boats are included.

2.2

ErgoClips

2.3 Do you ever weigh hazardous materials?

Your challenges

- Do you ever weigh hazardous materials? While handling samples, you might have to operate the balance in a glove box wearing gloves on both hands, as the material could be hazardous.
- Or you might weigh precious substances, where any loss of sample would directly affect costs.
- For weighing these special types of sample, the focus must be on sample handling for both safety and productivity.



SmartSens: Two optical sensors for hands-free operation are fitted to the XP balance terminal.

Our solutions

- Thanks to **"SmartSens"**, two optical sensors on the terminal of the **Excellence Plus XP**, basic weighing tasks can be performed without touching the balance.
- The Status bar shows which function of the sensor is activated by means of a green symbol on the terminal.
- For example: Zeroing, Taring the balance, Opening / Closing the draft shield (XP Analytical only), Printouts, M+ and IDs etc... can all be programmed to operate using the SmartSens.

We also offer the optional external sensor **"ErgoSens"**. It is possible to connect a maximum of two **ErgoSens** per Excellence balance. This means for the **XP balance**, a **maximum of 4 tasks** can be performed with just a wave of the hand. For **XS balances**, **maximum 2 tasks** can be performed at the same time (Tare, Zero or Printout).

Your benefits

- You can concentrate 100 % on sample handling.
- Contamination risk can be reduced to minimum.



ErgoSens: Optional sensor for hands-free operation can be connected using the Aux interface of the XP/XS balance.

Availability and ordering information

- **SmartSens** is mounted on the operating terminal of both **XP Analytical** and **Precision balances**.
- **ErgoSens** can be installed on **all Excellence balances** as an optional accessory (up to 2 ErgoSens per balance).
- A **Footswitch** for Aux interface is also available; the same operations as with the ErgoSens can be carried out by foot using the footswitch.

2.3

SmartSens & ErgoSens

Accessories	Description	Part No.
ErgoSens	Optional sensor for hands-free operation	11132601
Footswitch for Aux	Auxiliary footswitch with selectable function for the balance	11106741

2.4 Do you weigh in a hazardous environment?

Your challenges

- Your place of work is classified as an **Ex Zone 2** and you need a balance with **ATEX** approval to be able to perform weighings there.
- The size of your samples as well as the size of your containers fluctuates from small to large; the readability needed can vary from 0.1 mg to 1 g.
- You need to **print** the results of your formulations. However, you don't want to buy an ATEX printer and want to avoid too many cables.
- Your space is limited, your substance may be toxic, you would prefer to have the terminal of your balance separated from the platform.
- You need a **robust** balance with overload and IP protection.
- Your balance is used by various people: therefore it should be very easy to use.



Our solutions & your benefits

- We offer a complete range of **Excellence Precision balances for Ex Zone 2** (Category 3G) with 3 platform sizes (S, M and L) and a readability down to 0.1 mg.
- With the **Bluetooth** optional interface you can connect your balance wirelessly to a printer and/or a computer installed in a safe area, up to 10 meters from the balance (chapter 3.5).
- With the terminal on the **stand**, you do not only save space, you also avoid spilling substances onto the terminal.
- This robust balance offers complete **overload protection**. The soil-resistant design has **easy to clean** surfaces and edges.
- Thanks to the large **touch screen display**, operating the balance is easy and intuitive.
- Using an **ErgoSens** you can even operate the balance hands-free (chapter 2.3)
- For integration into existing systems, stand-alone platforms for Ex Zone 2 are also available.



Save space and avoid soiling the terminal using the optional stand.

► <http://www.mt.com/Ex>

Availability and ordering information

XS Ex-Balances

Models	Part No.
XS204SX	11121620
XS603SX	11121626
XS1003SX	11121632
XS5003SXDR	11121638
XS4002SX	11121644
XS6002SX	11121650
XS4001SX	11121656
XS6002MXDR	11121674
XS6001MX	11121680
XS12001MX	11121686
XS32001LX	11121704
XS64001LX	11121710
XS32000LX	11121716
XS64000LX	11121722

X Ex-Platform

Models	Part No.
X204SX	11121623
X603SX	11121629
X1003SX	11121635
X5003SXDR	11121641
X4002SX	11121647
X6002SX	11121653
X4001SX	11121659
X6001MX	11121683
X12001MX	11121689
X32001LX	11121707
X64001LX	11121713

It is possible to use the balance in hazardous areas while printing in safe areas without cable connections thanks to the BT optional interface.

Ex balances

2.5 Is your balance always adjusted correctly?

Your challenges

- Temperature changes may affect the sensitivity and sometimes the linearity of electronic balances, causing a difference between the displayed and true value.
- Periodical adjustments or tests at regular intervals may be required according to your SOP.


Two weights are installed with XP/XS Analytical balances, making **3 points linearity adjustment** possible (XP/XS Precision balances: 1 weight).

Our solutions for XP balances: proFACT & BalanceCheck

- Automatic adjustment using built-in reference weight.
- Temperature changes in the environment trigger the adjustment. Temperature criterion is selectable (0.5, 1, 2, 3 Kelvin).
- An adjustment record can be printed automatically, fulfilling GxP requirements.
- Automatic adjustment using a **built-in reference weight** and a prompt for Adjustment & Test using an **external weight** can be scheduled for a specific time (freely configurable: weekdays, up to 3 times per day).
- Up to 5 external weights for Adjustment & Test can be stored in the memory, and are easily selectable on the display.
- Adjustment & Test history events are recorded in the balance memory and can be printed (up to 50 events).

Our solutions for XS balances: FACT

- Automatic adjustment using **built-in reference weights** senses temperature changes in the environment.
- Adjustment record can be printed automatically fulfilling GxP requirements.
- Temperature criterion to trigger automatic adjustment is selectable.



proFACT & BalanceCheck: Periodical Adjustment & Test are implemented on schedule.

2.5

proFACT / FACT & BalanceCheck

Your benefits

- In the daily weighing process, your balance is always adjusted correctly. You are therefore freed from time-consuming maintenance duties.
- The Adjustment & Test history can be stored in the balance and printed anytime for traceability (proFACT).

Availability

- **proFACT** is included with **XP Analytical** and **Precision balances**
- **FACT** is included with **XS Analytical** and **Precision balances**

Function overview of Adjustment & Tests for XP and XS balances

	Automatic adjustment with internal weights		Prompting Adjustment & Test with external weights at a scheduled time	Number of programmable external weights for Adjustment & Test	Automatic print out of Adjustment & Test records	Adjustment & Test history record in the balance
	by temperature change	at a scheduled time				
XP	X	Weekdays, up to 3 times per day	Weekdays, up to 3 times per day	up to 5 weights	X	up to 50 records (the last 50)
XS	X	n/a	n/a	1 weight	X	n/a

2.6 Do you have cross-contamination risks?


Your challenges

- Powdery, liquid or gel samples are spilled on the draft shield or weighing pan.
- These samples could cause **cross-contamination**.
- Daily balance cleaning is time-consuming.

Our solutions for XP/XS Analytical balances

- The draft shield of the XP/XS Analytical balance is completely detachable.
- Spilled samples fall through the grid pan, and are retained on the drip tray below.
- Weighing pan: **SmartGrid** and drip tray can also be easily removed for cleaning.

Spilled samples are caught safely.



Completely detachable draft shield for the XP/XS analytical balance.



Every element can be detached.




The balance can be moved easily for benchtop cleaning.

2.6

**Cleaning
Concept for
Excellence
balances**

Our solutions for XP/XS Precision balances

- Straight-forward design for easy cleaning.
- A replaceable in-use protective cover comes as standard.
- **IP54 in use** for trouble-free weighing over many years.



Inner protection for spilled samples



Switches and connectors are sealed securely



Surfaces and edges designed for easy cleaning

Your benefits

- With our **Excellence** series, cleaning is easier than ever.
- The risk of cross-contamination can be minimized.
- **IP54 in use*** ensures longevity in validated processes and no downtime in production.
(* XP/XS Precision balances)



3. Increased Productivity

3

3.1 Would you like to increase the speed of your weighings?

Your challenges

- Air turbulence in the draft shield makes stabilization time slower (e.g. balance is in a laminar flow box).
- Both your hands are occupied with samples or tare containers. Automatic door opening and closing would enable more efficient weighing.

Our solutions & your benefits

- The innovative grid-shaped weighing pan, the "**SmartGrid**", prevents air turbulence from acting on the surface of the weighing pan, resulting in shorter and more stable weighing times. Repeatability is also improved.
- Optical sensor: **SmartSens** (on **XP Analytical balances**) opens the draft shield automatically with a simple wave of the hand.
- With **SmartSens** and the automatic door of the **XP Analytical**, there is no need to touch the balance to open the draft shield. Sampling speed is therefore faster than ever.
- Freely configurable door lever of the **XS Analytical** allows direct access of the sample onto the "SmartGrid" weighing pan.

Innovative **SmartGrid**



"Open sesame!" - automatic doors. Sampling speed is faster than ever (XP analytical balance).

Intelligent draft shield for XP analytical balance

Availability

- **SmartGrid** is fitted as standard on all **XP/XS Analytical balances**.
- **SmartSens** and **Automatic Doors** are fitted as standard on all **XP Analytical balances**.

Other advantages

- For faster weighing results, the less the draft shield is opened, the better the results. With **XP Analytical balances**, the amount that the door opens can be configured (depending on sample size; open 100 %, 75 %, 50 % or 25 %).
- For samples requiring a larger surface area on the weighing pan, a solid stainless steel plate is available for the **SmartGrid** (included as standard on all XP/XS Analytical balances).
- For toxic samples, disposable aluminum pans are also available for XP/XS Analytical balances. Part No.: 11106711 (includes 10 pcs).



Dominant hand is free: Right/Left door opening (XS analytical balance).



Disposable aluminium weighing pan.

3.1

**SmartGrid &
Draft Shield**

3.2 Would you like to improve the ease of use of the balance?

Your challenges

- It can be time-consuming to perform special weighing tasks (e.g. density measurement, statistical weighings, dynamic weighings etc...) with a conventional electronic balance.
- You would like to speed up the weighing process and to feel more confident using the balance.
- You have different balances used by people with different levels of experience. Therefore you need instruments with an intuitive and easy user interface.

Our solutions & your benefit

- The **Excellence balance range** offers simple touch screen operation and a dialogue oriented software guide enables you to perform even complicated applications easily and intuitively.
- Throughout the **Excellence balance range**, the user interface is the same. If you understand the operation of one Excellence balance, you can operate them all. **Know one, know all.** Less time for operator training.
- The major languages (English, French, German, Spanish, Italian, Japanese and Russian) are available in the software.
- For the **Excellence Plus XP series**, up to 8 display colors are selectable. You can identify which user settings are active at a glance.
- The XP's SmartTrac weighing-in aid has red and green indicators for fast and safe weighing.

Availability and functional overview


Series	Touch panel operation	7 languages	Color screen	Display size	Delta/SmartTrac
Excellence Plus XP	X	X	8 colors are selectable	117 x 87 mm	SmartTrac
Excellence XS	X	X	n/a	110 x 57 mm	DeltaTrac



Touch screen operation makes your weighing processes much easier.

3.2

Touch Screen Technology



The **touch screen** makes operation of the balance easy and safe.

3.3 Are you looking for software to manage weighing data?

Your challenges

- To manage and control a large amount of weighing data, you need powerful data management software for your laboratory.
- The software should be flexible, easy to use and it must ensure 100 % data security and traceability.

Our solutions & your benefits


With **LabX**, our laboratory instrument control software, you can perform the following tasks at the same time:

- **No mistakes in the data transfer:** Eliminate manual data entry and transcription errors.
- **Automatically transfer data to EXCEL:** Define the fields to be exported and the layout in LabX.
- **Seamless documentation:** Date and time stamp, multiple sample IDs, user-definable results and control limits, balance settings etc.
- **Secured and accessible data:** All the data is stored in an SQL database, so no data is lost.
- **Instrument management:** Test and adjustment scheduler, calibration history, weights library.
- **Integrated application:** Statistics, Differential weighing, Dynamic weighing, Interval weighing, Totalization etc...

For details about data security and LIMS integration, please refer to chapter 1.4.



LabX balance - data and instrument management software



Simple and flexible data transfer to EXCEL.

3.3

LabX




3.4 Do you weigh magnetic samples?

Your challenges

When magnetic samples are weighed, or if a magnetic stirring rod is weighed with samples, the weighing result can deviate from the true weight value.

Our solutions & your benefit

With a special **MPS** (**M**agnetic **P**rotection **S**hield) weighing pan, the magnetic effect can be reduced to a minimum. Even with magnetic samples, you can obtain reliable weighing results.



MPS (**M**agnetic **P**rotection **S**hield) weighing pan

Availability and ordering information

Accessories for the following balances	MPS weighing pan	Part No.
For XP / XS Precision S-platform models, with a readability of 1 g and 0.1 g	190 x 223 mm	11131542
For XP / XS Precision S-platform models, with a readability of 10 mg (also for DeltaRange models)	170 x 205 mm	11131541
For other models, individual measures required (e.g. below-the-balance weighing)	n/a	n/a

Remarks

- The MPS pan cannot be used on the XP1202S balance.
- The capacity of the balance is reduced slightly due to the heavier weighing pan.

3.4

MPS
Weighing
Pan

3.5 Are cables making a mess of your laboratory?

Your challenges

- Your work space is cluttered by too many cables and peripherals, or space is limited.
- You would like to connect your balance to a PC, printer or other peripherals at the same time.

Our solutions & your benefit

- With the optional **BTS** interface (**Bluetooth single point**), you can connect the balance to a PC, printer or other peripheral **wirelessly** up to 10 m. Cables are no longer necessary.
- With the optional **BT** interface (**Bluetooth multi point**), you can connect the balance to a PC, printer and other peripherals **wirelessly at the same time** (up to 6 peripherals).
- You can also isolate the balance in a laminar flow box or draft chamber, and acquire data by PC or printer wirelessly.

Availability and functional overview

Description		XP	XS	Part No.
Optional interface	BTS (single point)	X	X	11132535
Optional interface	BT (multi point)	X	X*	11132530
Wireless printer	BT-P42	X	X	11132540
Aux. wireless display	BT-BLD	X	X	11132555

* only available from terminal software version 3.00

MT's Bluetooth accessories are certified by the official Bluetooth qualification body.

For more information visit www.bluetooth.com

3.5

**Bluetooth
Wireless
Connection**




Bluetooth Interface

3.6 Do you need a customized draft shield?

For **XP/XS Precision balances**, we offer a variety of draft shields. If you need a larger or more functional draft shield for your XP/XS Precision balance, simply choose the best draft shield from the following table.

Availability overview

Description	Size (mm)	Optional availability						Part No.
								
		XP-S	XP/XS-S	XP/XS-M	XP/XS-L			
		0.1 mg	1 mg	10 mg	0.1 g			
ProDraft shield (standard for XP)	H: 247 mm	X	X	n/a	n/a	n/a	n/a	11131651
ProDraft shield, glass free (for food industry)	H: 247 mm	X	X	n/a	n/a	n/a	n/a	11131652
MagicCube draft shield (standard for XS)	H: 175 mm	n/a	X	n/a	n/a	n/a	n/a	11131650
Simple draft shield for 10 mg & 0.1 g models	H: 175 mm	n/a	n/a	X	X	n/a	n/a	11131653
Draft shield for whole balance (Medium)	W: 300 mm D: 450 mm H: 450 mm	n/a	X	X	X	X	n/a	11134430
Draft shield for whole balance (Large)	W: 550 mm D: 470 mm H: 580 mm	n/a	X	X	X	X	X	11134470

3.6

**Draft shields
for XP/XS
Precision
balances**

* To use Simple draft shield for XP/XS 0.1 g models, you need to order the weighing pan 170 x 205 mm (11132660).

3.7 Future software upgrades?

The software in our Excellence balances is continuously being developed.


Upgrades to the Excellence balance software can be performed easily and is free of charge using our "e-Loader" installation software.

e-Loader is available on our **Balance Support Website** and can be downloaded quickly and easily. Please visit and register on:

► <http://www.mt.com/balance-support-registration>

Ensures that your balance software is always the latest version!

For customers using the balance in a regulated environment, we recommend you ask your local MT technical support to install new software and help with its validation.



The online software update is easily accessible from your local PC with a standard internet connection.

3.7

e-Loader

Mettler-Toledo GmbH
Laboratory & Weighing Technologies
Address Im Langacher, CH-8606 Greifensee, Switzerland
Phone +41-1-944 36 36
Fax +41-1-944 31 90
E-mail LMSHotline@mt.com

METTLER TOLEDO

Text version 1.3

Supplier Declaration for Excellence Plus XP Software Update

Mettler-Toledo is certified according to ISO 9001. All products from Mettler-Toledo business unit Laboratory & Weighing Technologies are developed and manufactured based on the internationally recognized ISO 9001 quality system.

- The software is developed based on a life cycle model, using defined documentation, programming and test guidelines.
- All suppliers of system components must be formally assessed. Obviously standard products are an exception.
- The individual phases of the life cycle model and their working results are subject to a formal review (milestones).
- All documents and approved system components are subject to the modification check.

The platform software and terminal software for the XP balance line is in the Flash ROM on the motherboard on every balance. Both software types can be replaced if desired.

We confirm that the platform software version 4.11 and the terminal software 3.10 have been developed and tested according to the product development process described in the declaration of conformity. Important metrological data is stored separately in the EEPROM of the balance. This ensures that when the balance software is updated, metrological data cannot be modified intentionally or unintentionally.

The design is validated according to company-specific guidelines for the product and software development process. The relevant validation reports for these procedures and extensive documentation, such as source codes, can be viewed by arrangement.

On the support Website (<http://lab.mt.com/support>) you find two lists, one which describe the known and corrected errors up to Version 3.10 for the Terminal and up to Version 4.11 for the platform, and the other one that shows the new functionalities that are implemented in comparison to previous software. It should be noted that when new software versions are developed, only the modifications described in the Software Release Chart and their possible effects on the generation of the primary weight value are checked. We cannot accept liability for any secondary effects caused by software modification.

Mettler-Toledo GmbH
Laboratory & Weighing Technologies
Address Im Langacher, CH-8606 Greifensee, Switzerland
Phone +41-44-944 36 36
Fax +41-44-944 31 90
E-mail LMSHotline@mt.com

METTLER TOLEDO

Release date 1.12.2005

Software Release Chart - Excellence XS Analytical

Date	Application	Issue Description	Fixed in Terminal	Platform	TDNR* update needed
01.12.2005	Option Bluetooth Multipoint	Possibility to use a second interface Bluetooth multipoint (connection for up to 7 Bluetooth devices)	3.00	4.11	no
	Option MiniMettler	Possibility to use a second interface MiniMettler, for connection to older (legacy) systems			
	Option Ethernet	Possibility to use a second interface Ethernet for connection to Ethernet network			
	Option PS/2	Possibility to use a second interface PS/2 for connecting commercial keyboards and barcode readers			
	Piece counting	New application for piece counting.			
	% Weighing	New application for percent weighing.			
	Transfer key	This function key allows the current weight value to be sent directly to the connected host computer.			
	Languages	French, Spanish and Japanese language improved.			
	Statistic	Max n = 999 instead of 99			
	Display	Addition of special character.			
	Display	Graphic improved, when the screen save mode is switch on, the display stay inactive until the balance is ready.			
	Pretare with MT-SICS commands M06 and M26	When using the MT-SICS M06-command (Change of User/method number) and M26-command (Change of application) in conjunction with a PC application, the pretare value is transferred in error.			
	Protocol selection of ID	If an ID was chosen that ends with a dot (e.g. "Batch Nr."), this ID was not shown in the ID protocol selection menu nor in the info field			
1.10.2005	Stability of results	Filter optimization	1.14	4.11	no
11.06.2005	Stability of results	Enhanced robustness even under adverse environmental conditions	1.14	4.00	no
01.10.2004	Adjustment/Test protocol	On the adjustment protocol the date and time is always printed in the default format even if the format has been configured differently	1.14	1.13	no
	Display command (MT-SICS)	If a text is sent to the balance with the Display command that is too long to be shown on the balance screen then the text is only displayed partially. Instead of showing the least-significant text the most-significant text is shown.			
	Info Unit	If the "Custom Unit" is shown in the info field and later the "Custom Unit" is turned off, the "Custom Unit" will still be shown in the info field and in the set up menu the selection for "Info Unit" is empty.			
	Password protection	If the password protection for the system has been activated, the Date/Time can still be changed via touch display.			

* Only active with new TDNR (type definition number). To be updated by service technician only.

3.8 How soon can you expect to make a return on your investment?

You can appreciate that our **Excellence series** has the potential to enhance **your productivity and overall security** in weighing – security of the sample, of the result and of the user.

The following examples show a ROI calculation for the XS Analytical balance and XP Micro balance with optional accessories.

These examples should help you to estimate how fast your investment will pay for itself in your case. **Choosing an Excellence balance will certainly be a worthwhile investment!**



When speed is of the essence A supplier for the research-based chemical industry

Your challenges

- Produce single compounds for novel hits in HTS assays.
- Weigh a series of samples of a predefined volume according to customers' requirements.
- A lab technician fully dedicated to this task processes around 50 customer orders per day, with an average of 12 samples per order.
- The average weighing time contributes to 30% of the order. The speed of the balance is the key for higher throughput.

Our solution: Increase productivity through shorter weighing times

- **XS205DU:** Analytical balance with higher efficiency and productivity.
- **ErgoClip Basket** (Please refer to chapter 2.2)

Assumptions

- Compound orders produced/day: 50
- Total time/order: 9.6 min
- Weighing time/order: 3 min
- Margin/order sold: \$ 21.00

Return on investment: (using new XS205DU balance)

- 1 Productivity increase/day:
(20% higher weighing speed) $\frac{8 \text{ hours } (= 480\text{min})}{(3 \text{ min} \times 80\%) + 6.6 \text{ min}} - 50 = 3.3 \text{ orders more/day}$
- 2 Margin increase/day:
(annualized) $3 \text{ orders} \times \$ 21.00 \times 220 \text{ days} = \$ 13,860.00$
- 3 Initial investment:
(XS205DU) \$ 5,700.00*

* Price may vary according to the market.

3.8

ROI
Calculation
Example

Pay-back period: Less than 5 months

When less is more

A customer in the pharmaceutical industry

Your challenges

- Weigh reference substances for developing chemical reference solutions to calibrate HPLC systems. Some substances have a value of several thousand Euros per gram.
- Generate as many reference solutions as possible from these precious substances, according to USP regulation.
- To minimize the cost per reference solution, sample sizes must be as small as possible and equal to or greater than the minimal sample weight (according to USP) determined for the balance in use.
- Currently, a semi-micro balance with 100g capacity and 0.01mg readability is used. To eliminate transfer errors, small amounts of samples are transferred directly into bigger tare containers (tare weights are about 10 - 30g).

Our solution: Reduce the minimum sample weight

- **XP56 microbalance:** Micro balance with the highest capacity and security.
- **ErgoClip Flask micro** and **MinWeigh Door micro** (Please refer to 2.2).
- The **AntiStatic kit** for statically charged samples (Please refer to 2.1).
- **MinWeigh service** (please refer to 1.2).

Assumptions

- Minimum sample weight of the current analytical balance (Capacity 110g, Readability 0.01mg): 25mg
- Minimum sample weight of the XP56 microbalance (Capacity 52g, Readability 0.001mg): 8mg*
- Sample cost/mg (average): EUR 4.00
- Samplings/day (average): 2

* With smaller tare containers and optimal ambient conditions, the minimum sample weight value can be reduced down to 2.1mg, or even lower.

Return on investment: (using new XP56 balance)



- 1 Lower consumption of reference substance/day:
(8mg - 25mg = -17mg) $17\text{mg} \times 2 = 34\text{ mg}$
- 2 Cost savings:
(annualized) $34\text{mg} \times \text{EUR } 4.00 \times 220\text{ days} = \text{EUR } 29'920.00$
- 3 Initial investment:
(XP56 + Accessories & Service) $\text{EUR } 21'670.00^*$

* Price may vary according to the market.

Pay-back period: Less than 9 months

What if...? Excellence XP: "Your insurance for secure results".

Non-compliance with defined maximum limits can have disastrous consequences **in a regulated environment**. The following cases may occur in your laboratory and could cause significant losses.

The built-in warning functions of the XP range can prevent or minimize various risks, leaving you to concentrate on weighing with peace of mind.

What if...	Severe consequences	Our solution
The balance was inclined at some point, but measurements continued to be taken.	The SOP does not permit weighing without leveling; the data is regarded as "out of spec", and production has to be stopped...	LevelControl (Chapter 1.1)
Small samples were unintentionally measured below minimum weight .	Precious samples are lost, measurements must be repeated, but delivery of the sample material sometimes needs a couple of months.	MinWeigh (Chapter 1.2)
The weighing parameter settings, which can influence the precision of the weighing result, were accidentally changed by an operator.	Auditors found a disparity between the SOP and the actual settings on the balance. You are now asked to investigate the cause and effects...	User Management (Chapter 1.3)



4. More flexible Connectivity

4

4.1 Do you connect your balance with other peripherals?


On **all Excellence balances** the fixed RS232 interface as well as one of the interface options below can be used simultaneously.

7 different optional interfaces are available. Independent of this, the two outputs Aux1 and Aux2 can also be used.

Optional interfaces for Excellence series

Interface option	Part No.	Availability		Purpose
		XP	XS	
LocalCAN (LC)	11132505	X	X	• 2nd interface for connection of up to 5 LC instruments • For backwards compatibility (see chapter 4.2)
Bluetooth single point (BTS)*	11132535	X	X	• Bluetooth interface for wireless connection • Point to point connection, 1 balance with 1 peripheral
Bluetooth multi point (BT)*	11132530	X	X	• Bluetooth interface for wireless connection • Multi connection of 1 balance with up to 6 peripherals
RS232C (RS)	11132500	X	X	• 2nd interface for connection with RS232C peripherals • e.g.: One RS for Host, the other for a Printer or Titrator
PS/2	11132520	X	X	• 2 nd interface for a commercial keyboard or barcode reader
MiniMettler (MM)*	11132510	X	X	• 2 nd interface for backwards compatibility (see chapter 4.2)
Ethernet*	11132515	X	X	• 2 nd interface for LAN connection

* May not work with previous software versions. For further information please contact your local MT technical support.



4.1

Optional Interfaces for Excellence Series

Connectivity with peripherals

METTLER TOLEDO and third-party devices		Designation Part No.	LC	BTS	BT	RS	PS/2	MM	Ethernet	USB	Aux 1 & 2
	Barcode reader Information transmitted to the balance AC adapter 230 V (EU) AC adapter 115 V (USA)	BCR-RS 21900879 21900882 21900883	5			1					
		BCR-PS/2 21900880 BCR-PS/2Y 21900881					3				
	Printer GLP, GMP printer with setting possibilities for header with time, date, identification and serial number	BT-P42 11132540		2	2						
		RS-P42 11132555	5			1					
	Auxiliary display as benchtop version BLD = LCD display, backlight AD = Fluorescence display	BT-BLD 11132555		2	2						
	Auxiliary display as stand version BLDS = LCD display, backlight AD = Fluorescence display	RS/LC-BLDS 11132630	4			1					
	LC-IO Relay interface Input/output module for management of up to 8 dispensing or control systems (only for XP).	LC-IO 21202217	4			1					
	LV11 Automatic feeder Automatic feeder for small items	LV11 21900608	4			1					
	LC-FS Foot/hand switch For operation of the balance by hand or foot, 1 selectable function	LC-FS 229060	3								
	Footswitch for Aux 1 & 2 Auxiliary footswitch with selectable functions for the balance	Aux. Foot Switch 11106741									3
	ErgoSens (optical sensor) Auxiliary interface with selectable function for the balance. For hands-free operation.	ErgoSens 1132601									3
	Laptop		5	2	2	1			6	*	
	PC keyboard						3				
	Balances with MM interface							7			

*: RS/USB converter required (third-party device).

1: Cable 11101051 required

2: Wireless connection (Bluetooth)

3: Fixed cable

4: Cable 229116 (length = 5 m); additional cable: 229161 (length = 1 m), 229115 (length = 2 m)

5: It is possible to connect instruments having an RS interface with the LocalCan (LC) interface, using the cable LC-RS9: 229065 (length = 2 m)

6: LAN cable

7: Backwards compatibility to MiniMettler system.


4.1

**Connectivity
with
Peripherals**

MT-SICS (METTLER TOLEDO – Standard Interface Command Set)

We offer the **standardized bidirectional commands** for METTLER TOLEDO balance called **"MT-SICS"**. You can integrate your balance into a system, and can operate the balance remotely or automatically using bidirectional commands.

Basic commands (e.g. "Sending weighing data", "Tare the balance", "Zeroing" etc...) are identical among the entire METTLER TOLEDO balance ranges. For details of the bidirectional commands, please refer to the **"Operating instructions Excellence MT-SICS"** (visit www.mt.com and search using the keyword "MT-SICS").



Notes

For older generations of METTLER TOLEDO balances, certain bidirectional commands might not be supported by XP/XS balances, even with a LC or MM interface. To ensure compatibility, we recommend consulting your local MT technical support to verify that the commands which are used in your current system are fully supported in the **MT-SICS commands** of XP/XS balances.

LIMS integration and 21 CFR Part 11 compliance with LabX pro

LabX is the laboratory instrument control & data management software from METTLER TOLEDO for **Balances** and **Titration**. There are 3 levels of LabX: pro, light and direct. For more detailed information, please see chapter 1.4.



4.1

MT-SICS Commands

4.2 Are you about to replace your old balance?

Recommended replacement models

Your balance	Our recommendation
AE, AT, AX series	Excellence Plus XP Analytical series
AE, AG series	Excellence XS Analytical series
PM, PR, SR series	Excellence Plus XP Precision series
PM, PJ, PE, PG-S, SG series	Excellence XS Precision series

* Recommended models will vary depending on your specific weighing application or regulatory requirements. Please consult your local MT representative to choose the right balance for you.

Backward compatibility for system integration

When it comes to replacing your old METTLER TOLEDO balance which is integrated into your system, please make sure that the same type of interface is available. If your balance has a LocalCAN or a MiniMettler interface,

you may have to install optional interfaces in order to integrate the Excellence series into your existing system.


Your MT balance	Standard interface	Required optional interface for XP/XS	Part No.
AG, PG, PR, SR series ¹⁾ , etc...	LocalCAN (LC)	LC option for XP/XS	11132505
AT, PM, PJ, PE series ¹⁾ , etc...	MiniMettler (MM)	MM option for XP/XS ²⁾	11132510

1) The standard interface of each model may vary according to the production period. If it is difficult to identify the interface type of your balance, please ask your local MT technical support. 2) For the compatible functionality of MM option for XP/XS with older balance, please refer operating instructions or ask your local MT technical support.

Continuous mode

The **XP/XS** has 3 continuous modes* for **"MT-SICS"**, **"AT/MT"** and **"PM"**. You may use continuous mode with your AT or PM balance. If you replace an AT or PM with an XP/XS balance, you do not need to purchase the optional MiniMettler interface. In the peripherals setting menu, "Continuous mode" can be selected. The balance will send weighing data periodically and automatically, with an identical data output format than the AT or PM balance.

* Continuous mode: Sends weighing data periodically and automatically via Host mode.



On

MT-SICS

PM

AT/MT

3 different continuous modes are selectable

4.2

**Replacement
Models &
Backward
Compatibility**

4.3 Do you want to use a USB interface?

USB interface

Can computers with a **USB interface** be connected to METTLER TOLEDO **Excellence balances**?


Yes, all computers equipped with a USB interface can be connected to a METTLER TOLEDO Excellence balance using a "standard USB to RS converter", which is connected to the balance's RS232 interface. The full functionality of the classic RS-RS connection remains available.

Availability

USB to RS232 converters are country-specific products, requiring the installation of software on the computer. This software is supplied together with the converter. We recommend that the converter is purchased at your local electronic or PC supplier. We also recommend you consult your local MT technical support in order to check which converter has been tested locally and is therefore recommended for use with **Excellence balances**.

Example of connecting a laptop PC and a METTLER TOLEDO printer with a USB interface

The most common way to set up such a system is outlined below. Please refer to the drawing.



4.3

**USB
Interface**

Disclaimer

METTLER TOLEDO cannot guarantee the performance of third party products. It is therefore the responsibility of the customer to validate correct system operation and integration. Validation support and services are available from METTLER TOLEDO.

4.4 Do you connect your balance to a titrator?



Connection to titrator

A titrator is one of the most common instruments to be connected to a balance.

The following information shows how to configure **XP/XS balances** to be connected to METTLER TOLEDO Titrators.

1. Connection of XP/XS to actual MT-titrators

Titration	DL70ES/77		DL50/53/58		DL31/38 & DL32/39	
Software Version	3.1		2.4		2.12	
Titration Interface	CL20mA		RS232		RS232 fixed	
Communication	Unidirectional	Bidirectional	Unidirectional	Bidirectional	Unidirectional	Bidirectional
Configuration (fix)	2400, 7bit, even, no HS		2400, 7bit, even, no HS		2400, 7bit, even, no HS	

Interfaces used	Converter LC-CL5	Cable MM-CL5	Converter LC-CL5	Cable RS9-RS9	Converter LC-RS9	Cable RS9-RS9	Converter LC-RS9	Cable RS9-RS9	Converter LC-RS9	Cable RS9-RS9	Converter LC-RS9
Cable / Converter Part No.	ME229130		ME229130	ME11101051	ME229065	ME11101051	ME229065	ME11101051	ME229065		ME229065
Switch Settings	* 7 - 3 - 4 (continuous) (2400, 7/E, noHS)		* 0 - 3 - 4 (2400, 7/E, noHS)		* 7 - 3 - 4 (continuous) (2400, 7/E, noHS)		* 0 - 3 - 4 (2400, 7/E, noHS)		* 7 - 3 - 4 (continuous) (2400, 7/E, noHS)		* 0 - 3 - 4 (2400, 7/E, noHS)

Balance Interface	LC-Option	MM-Option	LC-Option	RS232 fix	LC-Option	RS232 fix or RS-Option	LC-Option	RS232 fix	LC-Option	RS232 fix or RS-Option	LC-Option
Configuration	Second display: LC-Option	Host: RS-Option	Host: LC-Option	Host: RS232 fix	Second display: LC-Option	Host: RS232 fix or RS-Option	Host: LC-Option	Host: RS232 fix	Second display: LC-Option	Host: RS232 fix or RS-Option	Host: LC-Option
Baudrate	* LC-Option	2400	* LC-Option	2400	* LC-Option	2400	* LC-Option	2400	* LC-Option	2400	* LC-Option
Parity		7bit/even/1		7bit/even/1		7bit/even/1		7bit/even/1		7bit/even/1	
Handshake		None		None		Xon/Xoff		None		Xon/Xoff	
Send mode		continuous		continuous		off		continuous		off	
Output-Format		AT		MT-SICS, PM, or AT				MT-SICS, PM, or AT			
Update/sec		5		5				5			
Weight obtained by	Cont.-mode	Cont.-mode	Titration	Cont.-mode	Cont.-mode	Titration	Titration	Cont.-mode	Cont.-mode	Titration	Titration

2. Connection of XS/XP to older MT-titrators

Titration	DL18/35	DL21/25	DL37
Software Version	1.2	1.2	3.2
Titration Interface	CL20mA	CL20mA	"Balance"
Communication	Unidirectional	Unidirectional	Unidirect. Bidirectional
Configuration (fix)	2400, 7bit, even, no HS	2400, 7bit, even, no HS	2400, 7bit, even, no HS

Interfaces used	CLC-CL5	Cable MM-CL5	LC-CL20	Cable MM-CL5	Cable Balance-RS	Converter LC-RS9
Cable / Converter Part No.	ME229130 + ME106090		ME229130 + ME106090		ME59759 + Gender Changer DB9 male-male	ME229065
Switch Settings	* 7 - 3 - 4 (continuous) (2400, 7/E, noHS)		* 7 - 3 - 4 (continuous) (2400, 7/E, noHS)			

Balance Interface	LC-Option	MM-Option	LC-Option	MM-Option	RS232 fix	RS232 fix or RS-Option	LC-Option
Configuration	Second display: LC-Option	Host: MM-Option	Second display: LC-Option	Host: MM-Option	Host: RS232 fix	Host: RS232 fix or RS-Option	Host: LC-Option
Baudrate	* LC-Option	2400	* LC-Option	2400	2400	2400	* LC-Option
Parity		7bit/even/1		7bit/even/1	7bit/even/1	7bit/even/1	
Handshake		None		None	None	Xon/Xoff	
Send mode		continuous		continuous	continuous		
Output-Format		AT		AT	AT		
Update/sec		5		5	5		
Weight obtained by	Cont.-mode	Cont.-mode	Cont.-mode	Cont.-mode	Cont.-mode	Titration	Titration

Notes: Weight transfer with Printkey (in Printer-mode) not recommended: Titration only accepts gross weights in grams. Other weight Units, Net-Weights or Weight results from approved balances are rejected by the titration.

Connection to a titration using LabX

With LabX, the instrument control software, you can integrate both **titrations** and **balances** into the same network, and can connect them to a LIMS for comprehensive data management. For more information, please refer to chapter 1.4.

* For setting up **XP/XS balances** with other titration brands, please contact your local MT technical support.

4.4

Titration Connection

[illegible]



5. Comprehensive Range of Accessories

5

5.1 Accessories overview for Excellence series

Description		Part No.	XP Analytical	XS	XP Precision	XS
Sensors	ErgoSens, optical Sensor for hands-free operation	11132601	X	X	X	X
Weighing pans	MPS (Magnetic Protection Shield) for 0.1 g models 190 mm x 223 mm	11132625	n/a	n/a	X	X
	MPS (Magnetic Protection Shield) for 10 mg models 170 mm x 205 mm	11132626	n/a	n/a	X	X
	Weighing pan 190 x 223 mm, incl. pan support	11132655	n/a	n/a	X	X
	Weighing pan 170 x 205 mm, incl. pan support and draft ring	11132660	n/a	n/a	X	X
Draft shields	Pro Draft Shield for 1 mg models, effective height 248 mm	11131651	n/a	n/a	X	X
	Draft Shield glass free for 0.1 mg and 1 mg models (food industry), effective height 248 mm	11131652	n/a	n/a	X	X
	Draft Shield for 0.1 g and 10 mg models, effective height 175 mm	11131653	n/a	n/a	X	X
	MagicCube draft shield for 1 mg models, effective height 175 mm	11131650	n/a	n/a	X	X
	Draft Shield for whole balance (S & M platform, dimensions (W x D x H) 300 x 450 x 450 mm	11134430	n/a	n/a	X	X
	Draft Shield for whole balance (L platform), dimensions (W x D x H) 550 x 470 x 580 mm	11134470	n/a	n/a	X	X
Printer	BT-P42 printer with wireless Bluetooth connection to balance	11132540	X	X	X	X
	RS-P42 printer with RS232C connection to balance	229265	X	X	X	X
Optional interfaces	BT option: 2nd interface Bluetooth, wireless multipoint connection for up to 6 Bluetooth devices	11132530	X	X	X	X
	BTS option: 2nd interface Bluetooth, wireless point-to- point connection	11132535	X	X	X	X
	Ethernet option: 2nd interface Ethernet for connection to Ethernet network	11132515	X	X	X	X
	E-link IP65 EB01: Ethernet connection to the e-link network with IP65 protection	11120003	X	X	X	X
	PS/2 option: 2nd interface for connecting commercial keyboards and barcode readers	11132520	X	X	X	X
	RS232C option: 2nd interface for connection of a printer (RS-P42), computer or titrator	11132500	X	X	X	X
	LocalCAN option: 2nd interface for connection of up to five LC (LocalCAN) instruments	11132505	X	X	X	X
	MiniMettler option: 2nd interface MiniMettler, for connection to older (legacy) systems	11132510	X	X	X	X

**XP/XS
Accessories**

5.1

Description		Part No.	XP Analytical	XS	XP Precision	XS
Protective covers	Protective cover for XP terminal (terminal only)	11132570	X	n/a	X	n/a
	Protective cover for XS terminal (terminal only)	11106870	n/a	X	n/a	X
	Protective cover for weighing platform "S", 10 mg / 0.1 g models (platform only)	11133034	n/a	n/a	X	X
	Protective cover for weighing platform "M" (platform only)	11132574	n/a	n/a	X	X
	Protective cover for XS balance with "S" platform, 10 mg / 0.1 g models (whole balance)	11132571	n/a	n/a	n/a	X
	Protective cover for XS balance with "M" platform (whole balance)	11132572	n/a	n/a	n/a	X
	Protective cover for XS terminal attached to a "L" platform (terminal only)	11132573	n/a	n/a	n/a	X
Stands and wall mounting	Terminal stand for placement of the terminal 30 cm above weighing pan (S & M Platform)	11132636	n/a	n/a	X	X
	Terminal stand for placement of the terminal 30 cm above weighing pan (L Platform)	11132653	n/a	n/a	X	X
	Stand for Terminal and Printer	11106730	X	X	n/a	n/a
	Terminal wall mounting	11132665	X	X	X	X
ErgoClips	ErgoClip Basket small	11106747	X	X	n/a	n/a
	ErgoClip Titration Basket	11106883	X	X	n/a	n/a
	ErgoClip Weighing Boat	11106748	X	X	n/a	n/a
	ErgoClip Round Bottom Flask	11106746	X	X	n/a	n/a
	ErgoClip Flask	11106764	X	X	n/a	n/a
	MinWeigh Door (ideal to use with "ErgoClip Flask")	11106749	X	X	n/a	n/a
	ErgoClip Flask micro (for XP micro balances only)	11107879	X	n/a	n/a	n/a
	MinWeigh Door micro (for XP micro balances only, ideal to use together with "ErgoClip Flask micro")	11107869	X	n/a	n/a	n/a
	ErgoClip Tube	11106784	X	X	n/a	n/a
	Single-use aluminum weighing pan (10 pans)	11106711	X	X	n/a	n/a
	Weighing Kit - Set of ErgoClips (Basket small, Weighing Boat, Round Bottom Flask and Single use aluminum pans are included)	11106707	X	X	n/a	n/a
AC adapter protection	IP54 AC adapter protection cover	11132550	X	X	X	X
Anti-theft device	Steel cord	11600361	X	X	X	X
Terminal cable	Terminal extension cable, 4.5 m	11600517	X	X	X	X
Density kit	Density kit for precision balances, 1 mg and 0.1 mg models	11132680	n/a	n/a	X	X
	Density kit for all analytical balances	11106706	X	X	n/a	n/a
	Certified thermometer for density determination	11132685	X	X	X	X
	Sinker 10 ml: for density determination of liquids	210260	X	X	X	X
	Sinker 10 ml, certified model	210672	X	X	X	X
Dynamic weighing kit	Dynamic weighing kit for 0.1 mg and 10 mg models, 4 liter bowl	11132657	n/a	n/a	X	X

XP/XS Accessories

5.1

Description		Part No.	XP Analytical	XS	XP Precision	XS
Auxiliary displays	BT-BLD Bluetooth aux. display for benchtop mounting. Wireless connection. Backlit LCD display.	11132555	X	X	X	X
	RS/LC-BLDS auxiliary display for benchtop mounting or on balance. Backlit LCD display.	11132630	X	X	X	X
	LC-AD auxiliary display on benchtop stand with fluorescent display	229140	X	X	X	X
	LC-ADS auxiliary display on high benchtop stand with fluorescent display	229150	X	X	X	X
	LC/RS-BLD auxiliary display on benchtop stand with backlit LC display	224200	X	X	X	X
Barcode reader	RS232C barcode reader	21900879	X	X	X	X
	AC adapter for barcode reader 230 V EUR	21900882	X	X	X	X
	AC adapter for barcode reader 115 V USA	21900883	X	X	X	X
	PS/2 barcode reader	21900880	X	X	X	X
	PS/2Y barcode reader	21900881	X	X	X	X
Filling-process control	LV11 small-item transporter for automatic loading of small items on the balance	21900608	X	X	X	X
	LV11 Draft shield door "Pro" for 0.1 mg and 1 mg models	11132711	n/a	n/a	X	X
	LV11 Draft shield door for all analytical balances	11106715	X	X	n/a	n/a
LC-I/O Box	Relay interface for control of up to 8 external instruments from the balance	21202217	X	n/a	X	n/a
LC-Footswitch	Auxiliary footswitch with selectable function for the balance	11106741	X	X	X	X
	LC-FS with selectable function	229060	X	X	X	X
	Switchbox, connect up to 3 balances with LocalCan interface to a printer	229220	X	X	X	X
Transport case	Transport case for precision balances, 10 mg / 0.1 g models. With space for printer.	11132595	n/a	n/a	X	X
	Transport case for all analytical balances	11106873	X	X	n/a	n/a
XP-SE-kit	Separate electronic kit for XP analytical balances	11106743	X	n/a	n/a	n/a
	Extension cable 0.6 m	211535	X	n/a	n/a	n/a
	Extension cable 5 m	210688	X	n/a	n/a	n/a
PC Volume	PC Volume software (Option 1)	21900791	X	X	n/a	n/a
	Evaporation trap and PC Volume software (Option 3)	21901263	X	X	n/a	n/a
	Evaporation trap (Option 2)	21901264	X	X	n/a	n/a
AntiStatic kit	Built-in AntiStatic Kit for XP analytical balances, complete with power supply	11107761	X	n/a	n/a	n/a
	Optional second point-electrode for AntiStatic Kit	11107762	X	n/a	n/a	n/a
	Universal AntiStatic kit complete (U-kit)	11107767	X	X	X	X
	Optional U-electrode for universal AntiStatic kit	11107764	X	X	X	X
	Optional point-electrode for universal AntiStatic kit	11107765	X	X	X	X
	Power supply for universal AntiStatic kit (requires optional electrode)	11107766	X	X	X	X
Production certificate	Pro production certificate for XS balances	11106895	incl.	X	incl.	X

XP/XS Accessories

5.1

Description			Part No.	XP Analytical	XS	XP Precision	XS
Hook below the balance	Hook for weighing below balance with "M" or "L" platform (integral to "S" platform)		11132565	n/a	n/a	X	X
LabX balance	LabX pro balance	English	11120301	X	X	X	X
		German	11120302	X	X	X	X
		French	11120303	X	X	X	X
	LabX light balance	English	11120317	X	X	X	X
		German	11120318	X	X	X	X
		French	11120319	X	X	X	X
	LabX direct balance	English	11120340	X	X	X	X
	Validation handbook I	English	11120332	X	X	X	X
	Validation handbook II	English	11120333	X	X	X	X

**XP/XS
Accessories**

5.1

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and extend across the width of the page. There are no margins, text, or other markings on the paper.

Mettler-Toledo GmbH,
Laboratory & Weighing Technologies
 CH-8606 Greifensee, Switzerland
 Phone +41-44-944 22 11
 Fax +41-44-944 30 60

Australia	Mettler-Toledo Ltd., Port Melbourne, Victoria 3207 Tel. (03) 9644 57 00, Fax (03) 9645 39 35
Austria	Mettler-Toledo GmbH., A-1100 Wien Tel. (01) 604 19 80, Fax (01) 604 28 80
Belgium	Mettler-Toledo s.a., B-1932 Zaventem Tel. (02) 334 02 11, Fax (02) 334 03 34
Brazil	Mettler-Toledo Ltda., 06455-000 Barueri/São Paulo Tel. (11) 7295 1692, Fax (11) 421 3459
Canada	Mettler-Toledo Inc., Ontario, Canada Tel (800) 638-8537, Fax (905) 681-8036
China	Mettler-Toledo (Shanghai) Ltd., Shanghai 200233 Tel. (21) 6485 04 35, Fax (21) 6485 33 51
Croatia	Mettler-Toledo, d.o.o., HR-10000 Zagreb Tel. 01 29 58 130, Fax 01 29 58 140
Czech Republic	Mettler-Toledo, spol. s.r.o., CZ-100 00 Praha 10 Tel. 272 123 150, Fax 272 123 170
Denmark	Mettler-Toledo A/S, DK-2600 Glostrup Tel. (43) 270 800, Fax (43) 270 828
France	Mettler-Toledo s.a., F-78222 Viroflay Tel. (01) 3097 17 17, Fax (01) 3097 16 16
Germany	Mettler-Toledo GmbH, D-35353 Giessen Tel. (0641) 50 70, Fax (0641) 507 128
HongKong	Mettler-Toledo (HK) Ltd., Kowloon Tel. (852) 2744 1221, Fax (852) 2744 6878
Hungary	Mettler-Toledo, KFT, H-1139 Budapest Tel. (01) 288 40 40, Fax (01) 288 40 50
India	Mettler-Toledo India Pvt Ltd, Mumbai 400 072 Tel. (22) 2857 08 08, Fax (22) 2857 50 71
Italy	Mettler-Toledo S.p.A., I-20026 Novate Milanese Tel. (02) 333 321, Fax (02) 356 29 73
Japan	Mettler-Toledo K.K., Tokyo 143 Tel. (3) 5762 0606, Fax (3) 5762 0756
Korea	Mettler-Toledo (Korea) Ltd., Seoul 137-130 Tel. 82-2-3498-3500, Fax 82-2-3498-3555
Kazakhstan	Mettler-Toledo Central Asia, 480009 Almaty Tel. (07) 3272 980 834, Fax (07) 3272 980 835
Malaysia	Mettler-Toledo (M) Sdn. Bhd., 40150 Selangor Tel. (603) 784 55 773, Fax (603) 784 58 773
Mexico	Mettler-Toledo S.A. de C.V., México C.P. 06430 Tel. (55) 5547 5700, Fax (55) 5547 2128
Netherlands	Mettler-Toledo B.V., NL-4004 JK Tiel Tel. (0344) 638 363, Fax (0344) 638 390
Norway	Mettler-Toledo A/S, N-1008 Oslo Tel. 22 30 44 90, Fax 22 32 70 02
Poland	Mettler-Toledo, Sp. z o.o., PL-02-822 Warszawa Tel. (22) 545 06 80, Fax (22) 545 06 88
Russia	Mettler-Toledo C.I.S. AG, 10 1000 Moskau Tel. (095) 921 92 11, Fax (095) 921 78 68
Singapore	Mettler-Toledo (S) Pte. Ltd., Singapore 139959 Tel. 65-6890 0011, Fax 65-6890 0012
Slovakia	Mettler-Toledo, SK-82104 Bratislava Tel. (02) 434 27 496, Fax (02) 433 37 190
Slovenia	Mettler-Toledo, d.o.o., SI-1236 Trzin Tel. (01) 530 8050, Fax (01) 562 1789
Spain	Mettler-Toledo S.A.E., E-08908 Barcelona Tel. (93) 223 76 00, Fax (93) 223 76 01
Sweden	Mettler-Toledo AB, S-12008 Stockholm Tel. (08) 702 50 00, Fax (08) 642 45 62
Switzerland	Mettler-Toledo (Schweiz) AG, CH-8606 Greifensee Tel. (044) 944 45 45, Fax (044) 944 45 10
Thailand	Mettler-Toledo (Thailand) Ltd., Bangkok 10320 Tel. 66-2-723 0300, Fax 66-2-719 6479
Taiwan	Mettler-Toledo Pac Rim AG, Taipei Tel. (2) 2657 8898, Fax (2) 2657 0776
United Kingdom	Mettler-Toledo Ltd., Leicester, LE4 1AW Tel. (0116) 235 70 70, Fax (0116) 236 5500
United States	Mettler-Toledo, Inc., Columbus, OH 43240 Tel. 800-METTLER, Fax (614) 438 4900

For all other countries

Mettler-Toledo GmbH
 PO Box VI-400, CH-8606 Greifensee
 Tel. +41-44-944 22 11, Fax +41-44-944 31 70

Subject to technical changes
 © 03/2006 Mettler-Toledo GmbH
 Printed in Switzerland 11795638

www.mt.com/excellence

More information



Quality certificates

Development, production and testing under ISO 9001.



Environment management system under ISO 14001.



"Communauté Européenne".

This mark assures you that our products comply with the latest guidelines.